

Anti-SPO11 antibody (269-358 aa) (STJ11107536)

STJ11107536

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	23626
Gene Symbol	SPO11
Uniprot ID	SPO11_HUMAN
Immunogen	
Immunogen Region	269-358 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 269-358 of human SPO11 (NP_937998.1).
Immunogen Sequence	AIRWLGLLPDDLKRLNVPKD SLIPLTKRDQMKLDSILRRP YVTCQPFWRKEMEIMADSKM KAEIQALTFLLSSDYLSRVYL PNKLKFGGWI

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081