

Anti-CACNA1C antibody (1780-1980 aa) (STJ11107538)

STJ11107538

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	775
Gene Symbol	CACNA1C
Uniprot ID	CAC1C_HUMAN
Immunogen	
Immunogen Region	1780-1980 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1780-1980 of human CACNA1C (NP_955630.3).
Immunogen Sequence	KAGSSQGDTEPSHEKLVDS TFTPSSYSSTGSGNANINNAN NTALGRLPRPAGYPSTVSTV EGHGPPLSPAIRVQEVAVKL SSNRERHVPKCEDLELRDS GSAGTQAHCLLLRRANPSRC HSRESQAAMAGQEETSQDET YEVKMNHDTEACSEPSLLST EMLSYQDDENRQLTLPEEDK RDIRQSPKRGFLRSASLGRR A