

Anti-RNF135 antibody (121-309 aa) (STJ11107541)

STJ11107541

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:1000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	84282
Gene Symbol	RNF135
Uniprot ID	RN135_HUMAN
Immunogen	
Immunogen Region	121-309 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 121-309 of human RNF135 (NP_115698.3). ELQRRVAVEKSITEVAQELTE LVEHLVDIVRSLQNQRPLSE SGPDNELSLGKAFSSGVDL SMASPKLVTSDDTAAGKIRDI LHDLEEIQEKLQESVTWKEA PEAQMQGELLEAPSSSSCPL PDQSHPALRRASRFAQWAIH PTFNLKSLSCSLEVSKDSRT VTVSHRPQPYRWSCERFSTS QVLCSSQALS
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081