

Anti-IDUA antibody (354-653 aa) (STJ11107544)

STJ11107544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3425
Gene Symbol	IDUA
Uniprot ID	IDUA_HUMAN
Immunogen	
Immunogen Region	354-653 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 354-653 of human IDUA (NP_000194.2).
Immunogen Sequence	SYHPPHFAQRTLTARFQVNN TRPPHVQLLRKPVLTAMGLL ALLDEEQLWAEVSQAGTVLD SNHTVGVLASAHRPQGPADA WRAAVLIYASDDTRAHPNRS VAVTLRLRGVPPGPGLVVYT RYLDNGLCSPDGEWRRRLGRP VFPTAEQFRMRRAAEDPVAA APRPLPAGGRLTLRPALRLP SLLLVHVCARPEKPPGQVTR LRALPLTQGQLVLVWSDEHV GSKCLWYEQFSQDGKAY