

Anti-eIF4EBP3 antibody [SMM] (STJ11107551)

STJ11107551

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Transfected

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG1k
Storage Instruction	

TARGET INFORMATION

Gene ID	8637
Gene Symbol	EIF4EBP3
Uniprot ID	4EBP3_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human eIF4EBP3
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081