

Anti-GST-Tag antibody (1-218 aa) [SMM] (STJ11107564)
STJ11107564

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IP/ELISA
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2000-1:5000
Range	IP:3 Mu g antibody for 1 Mu g extracts of recombinant protein ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 1.5% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG1
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	1-218 aa
Specificity	Recombinant protein containing a sequence corresponding to amino acids 1-218 of GST protein.
Immunogen Sequence	MSPILGYWKIKGLVQPTRLL LEYLEEKYEELHYERDEGDK WRNKKFELGLEFPNLPYYID GDVKLTQSMAIRYIADKHN MLGGCPKERAIEISMLEGAVL DIRYGVSR IAYSKDFETLKV DFLSKLPEMLKMFEDRLCHK TYLNGDHVTHPDFMLYDALD VLYMDPMCLDAFPKLVCFK KRIEAIPIQIDKYLKSSKYIA WPLQGWQATFGGGDHPPK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081