

Anti-His-Tag antibody [SMM] (STJ11107566)

STJ11107566

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/IP/ELISA
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:20000 IF/ICC:1:50-1:200 IP:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with pH 7.3.
Isotype	IgG1
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	HHHHHHH
Immunogen Region	
Specificity	A synthetic peptide corresponding to His tag.
Immunogen Sequence	HHHHHHH