

Anti-GST-Tag antibody [SRM] (STJ11107590)

STJ11107590

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IP
Host/Source	Rabbit
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IP:1:500-1:5000
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human GST-Tag.
Immunogen Sequence	