

Anti-MYC-Tag antibody [SRM] (STJ11107591)
STJ11107591**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	IP
Host/Source	Rabbit
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IP:1:50-1:200
Formulation	
Isotype	
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	A synthetic peptide corresponding to Myc-Tag.
Immunogen Sequence	