

Anti-GFP-Tag antibody (1-238 aa) [SRM] {Affinity Gel} (STJ11107594)

STJ11107594

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IP
Host/Source	Rabbit
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Affinity Gel
Purification	Affinity purification
Dilution Range	IP:20 Mu l-40 Mu l Affinity Gel for 100 Mu g-300 Mu g extracts of whole cells
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	1-238 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of GFP.
Immunogen Sequence	MSKGEELFTGVVPIVELDG DVNGHKFSVSGEGEGDATYG KLTCLKICTTGKLPVPWPTL VTTLTYGVCFSRYPDHMKR HDFFKSAMPEGYVQERTISF KDDGNYKTRAEVKFEGDTLV NRIELKGIDFKEDGNILGHK LEYNYNSHNVYITADKQKNG IKANFKTRHNIEDGGVQLAD HYQNTPIGDPVLLPDNHY LSTQSALSKDPNEKRDHMLV LEFVTAAGITHGMDLYK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081