

Anti-Phospho-IKK Beta-Y199 antibody (STJ11107629)

STJ11107629

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3551
Gene Symbol	IKBKB
Uniprot ID	IKKB_HUMAN
Immunogen	
Immunogen Region	
Specificity	A phospho specific peptide corresponding to residues surrounding Y199 of human IKK Beta
Immunogen Sequence	