

Anti-Phospho-Tyrosine Hydroxylase-S19 antibody (STJ11107634)

STJ11107634

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7054
Gene Symbol	TH
Uniprot ID	TY3H_HUMAN
Immunogen	
Immunogen Region	
Specificity	A phospho specific peptide corresponding to residues surrounding S19 of human Tyrosine Hydroxylase
Immunogen Sequence	