

Anti-Phospho-IGF1R-Y1135/Y1136 + Y1150/Y1151 antibody (STJ11107647)

STJ11107647

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3643 3480
Gene Symbol	INSR IGF1R
Uniprot ID	INSR_HUMAN IGF1R_HUMAN
Immunogen	TDYYRK
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y1135 & Y1136 of human IGF1R/IGF1R (NP_000866.1).
Immunogen Sequence	TDYYRK