

Anti-Phospho-p90Rsk/RSK1/RPS6KA1-T573 antibody (STJ11107649)

STJ11107649

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6195
Gene Symbol	RPS6KA1
Uniprot ID	KS6A1_HUMAN
Immunogen	LMTPC
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around T573 of human p90Rsk/RSK1/RPS6KA1 (NP_002944.2).
Immunogen Sequence	LMTPC