

Anti-Phospho-MUC1-Y1203 antibody (STJ11107652)

STJ11107652

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4582
Gene Symbol	MUC1
Uniprot ID	MUC1_HUMAN
Immunogen	DTYHP
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y1203 of human MUC1 (P15941).
Immunogen Sequence	DTYHP