

Anti-Phospho-PI3KP85 Alpha/P55 Gamma/P85 Beta-Y467/Y199/Y464 antibody (STJ11107668)

STJ11107668

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5295 8503 5296 PIK3R1 PIK3R3 P85A_HUMAN P55G_HUMAN P85B_HUMAN
Immunogen	DRLYEE
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y467/Y199/Y464 of human p85-alpha (NP_852664.1).
Immunogen Sequence	DRLYEE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081