

Anti-Phospho-p53-S20 antibody (STJ11107670)

STJ11107670

GENERAL INFORMATION

| | |
|--------------------------|--------------------|
| Product Type | Primary antibodies |
| Short Description | |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | Lot specific |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3. |
| Isotype | IgG |
| Storage Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 7157 |
| Gene Symbol | TP53 |
| Uniprot ID | P53_HUMAN |
| Immunogen | TFSDL |
| Immunogen Region | |
| Specificity | A synthetic phosphorylated peptide around S20 of human p53 (NP_000537.3). |
| Immunogen Sequence | TFSDL |