

Anti-Phospho-Rptor-S791 antibody (STJ11107673)

STJ11107673

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	Rptor
Uniprot ID	RPTOR_MOUSE
Immunogen	
Immunogen Region	
Specificity	A phospho specific peptide corresponding to residues surrounding S791 of mouse Rptor
Immunogen Sequence	