

Anti-Phospho-DOK7-Y406 antibody (STJ11107689)

STJ11107689

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	231134
Gene Symbol	Dok7
Uniprot ID	DOK7_MOUSE
Immunogen	
Immunogen Region	
Specificity	A phospho specific peptide corresponding to residues surrounding Y406 of Mouse DOK7
Immunogen Sequence	