

Anti-Phospho-RNF168-S60 antibody (STJ11107691)

STJ11107691

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	165918
Gene Symbol	RNF168
Uniprot ID	RN168_HUMAN
Immunogen	VSSWT
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around S60 of human RNF168 (NP_689830.2).
Immunogen Sequence	VSSWT