

Anti-Phospho-RIG antibody (STJ11107693)

STJ11107693

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	230073
Gene Symbol	Rigi
Uniprot ID	RIGI_MOUSE
Immunogen	VRYKT
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y880 of mouse RIG-1 (NP_766277.3).
Immunogen Sequence	VRYKT