

Anti-Phospho-DDR1-Y513 antibody (STJ11107701)

STJ11107701

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	780
Gene Symbol	DDR1
Uniprot ID	DDR1_HUMAN
Immunogen	PAYRL
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y513 of human DDR1 (NP_054700.2).
Immunogen Sequence	PAYRL