

Anti-Phospho-NBS1/NBN-S343 antibody (STJ11107709)

STJ11107709

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4683
Gene Symbol	NBN
Uniprot ID	NBN_HUMAN
Immunogen	SLSQG
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around S343 of human NBS1/NBN (NP_002476.2).
Immunogen Sequence	SLSQG