

Anti-Phospho-PI3 Kinase p85 alpha-Y607 antibody (STJ11107710)

STJ11107710

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5295
Gene Symbol	PIK3R1
Uniprot ID	P85A_HUMAN
Immunogen	DQYSL
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y607 of human PI3 Kinase p85 alpha (NP_852664.1).
Immunogen Sequence	DQYSL