

Anti-Phospho-RIPK1/RIP-S166 antibody (STJ11107720)

STJ11107720

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	19766
Gene Symbol	Ripk1
Uniprot ID	RIPK1_MOUSE
Immunogen	TWSKL
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around S166 of mouse Phospho-RIPK1-S166 (NP_033094.3).
Immunogen Sequence	TWSKL