

Anti-Phospho-BRD4-S1083 antibody (STJ11107721)

STJ11107721

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	23476
Gene Symbol	BRD4
Uniprot ID	BRD4_HUMAN
Immunogen	ASTPP
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around S1083 of human Phospho-BRD4-S1083 (NP_490597.1).
Immunogen Sequence	ASTPP