

## Anti-SCFD1 antibody (373-642) (STJ111440)

STJ111440

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IHC-P/ELISA     |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                            |                                                                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                                                 |
| <b>Clone ID</b>            |                                                                                                                                                            |
| <b>Concentration</b>       | Lot specific                                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                               |
| <b>Purification</b>        | Affinity purification                                                                                                                                      |
| <b>Dilution</b>            | WB:1:500-1:2000                                                                                                                                            |
| <b>Range</b>               | IHC-P:1:50-1:200<br>ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                                                                        |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                             |

### TARGET INFORMATION


|                           |                                                                                                                                                                                                                                                              |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 23256                                                                                                                                                                                                                                                        |
| <b>Gene Symbol</b>        | SCFD1                                                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>         | SCFD1_HUMAN                                                                                                                                                                                                                                                  |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                              |
| <b>Immunogen Region</b>   | 373-642                                                                                                                                                                                                                                                      |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 373-642 of human SCFD1 (NP_057190.2).                                                                                                                                          |
| <b>Immunogen Sequence</b> | SMLSDNTAKLTSAVSSLPEL LEKKRLIDLHTNVATAVLEH IKARKLDVYFEYEEKIMSKT TLDKSLDIISDPDAGTPED<br>KMRLFLIYYISTQQAPSEAD LEQYKKALTDAGCNLNPLQY IKQWKAFKMASAPASYGST TTKPMGLLSRVMTGSGQFVM<br>EGVKNLVKQQNLVPTRILD NLMEMKSNPETDDYRYFDPK MLRGNDSSVPRNKNPFQEI VFVVGNGNYIEYQNLVDYI |




Western blot analysis of extracts of various cell lines, using SCFD1 antibody (STJ111440) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 $\mu$ g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry analysis of paraffin-embedded rat lung using SCFD1 Rabbit polyclonal antibody (STJ111440) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry of paraffin-embedded human liver cancer using SCFD1 rabbit polyclonal antibody (STJ111440) at dilution of 1:200 (40x lens).



Immunohistochemistry analysis of paraffin-embedded mouse kidney using SCFD1 Rabbit polyclonal antibody (STJ111440) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081