

Anti-NISCH antibody (500-700) (STJ111485)

STJ111485

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11188
Gene Symbol	NISCH
Uniprot ID	NISCH_HUMAN
Immunogen	
Immunogen Region	500-700
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human NISCH (NP_009115.2).
Immunogen Sequence	SLPQPILSNQGIMFVQEEAL ASSLSSTDLSLTPEHQPIAQG CSDSLESIPAGQAASDDL RD VPGAVGGASPEHAEPEVQVV PGSGQIIIFLPFTCIGYTATN QDFIQLRLTLIRQAIERQLP AWIEAANQREEGQGEGQEEEE DEEEEEEDVAENRYFEMGP PDVEEEEGGGQEEEEEEEE DEEAEEERLAEWALGADED F