

Anti-SLCO1C1 antibody (450-554) (STJ111626)

STJ111626

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	53919
Gene Symbol	SLCO1C1
Uniprot ID	SO1C1_HUMAN
Immunogen	
Immunogen Region	450-554
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 450-554 of human SLCO1C1 (Q9NYB5).
Immunogen Sequence	SDVAGLTVSYQGTPVSYHE RALFSDCNSRCKCSETKWEP MCGENGITYVSACLAGCQTS NRSQKNIIFYNCTCVGIAAS
Sequence	KSGNSSGIVGRCQKDNGCPQ MFLYF