

Anti-RGS9 antibody (230-350) (STJ111658)

STJ111658

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8787
Gene Symbol	RGS9
Uniprot ID	RGS9_HUMAN
Immunogen	
Immunogen Region	230-350
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 230-350 of human RGS9 (NP_003826.2).
Immunogen Sequence	MYQQALMRSTVKSSVSLGG IVKYSEQFSSNDMSGCLP SNPWITDDTQFWDLNAKLVE IPTKMRVERWAFNFSLEIRD PKGRQSFQYFLKKEFSGENL GFWEACEDLKYGDQSKVKEK A