

Anti-STK25 antibody (297-426) (STJ111804)

STJ111804

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10494
Gene Symbol	STK25
Uniprot ID	STK25_HUMAN
Immunogen	
Immunogen Region	297-426
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 297-426 of human STK25 (NP_006365.2).
Immunogen Sequence	EESSSEDSIDGAEDEGEQ PIWTFPPTIRPSPHSKLHKG TALHSSQKPAEPVKRQPRSQ CLSTLVRPVFGELKEKHKQS GGSGVGALEELNAFLAEES CPGISDKLMVHLVERVQRFS HNRNHLTSTR



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081