

Anti-NFYC antibody (27-335) (STJ111874)

STJ111874

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4802
Gene Symbol	NFYC
Uniprot ID	NFYC_HUMAN
Immunogen	
Immunogen Region	27-335
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 27-120 of human NFYC (NP_055038.2).
Immunogen Sequence	MEEIRNLTVKDFRVQELPLA RIKKIMKLEDEVKMISAEAP VLFKAAQIFITELTLRAWI HTEDNKRRTLQRNDIAMAIT KFDQDFDLIDIVPR



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081