

Anti-NRIP1 antibody (939-1158) (STJ111996)
STJ111996

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8204
Gene Symbol	NRIP1
Uniprot ID	NRIP1_HUMAN
Immunogen	
Immunogen Region	939-1158
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 939-1158 of human NRIP1 (NP_003480.2). LLLSENCVRDLSPHRNSVA DSKKKGHKNNVTNSKPEFSI SSLNGLMYSSTQPSSCMDNR TFSYPGVVKTVPVPTFFPEHL
Immunogen Sequence	GCAGSRPESGLLNGCSPMPSE KGPIKWITDAEKNEYEKDS PRLTKTNPILYYMLQKGGNS VTSRETQDKDIWREASSAES VSQVTAKEELLPTAETKASF FNLRSPYNHMGNNASRPHS ANGEVYGLLGSVLTIKKESE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081