

Anti-CSPP1 antibody (962-1221) (STJ112024)
STJ112024

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	79848
Gene Symbol	CSPP1
Uniprot ID	CSPP1_HUMAN
Immunogen	
Immunogen Region	962-1221
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 962-1221 of human CSPP1 (NP_079066.5).
Immunogen Sequence	QSPKGLDAATFQNVHDFNEL KDRDSETRVDLKFMYLDPPR DHHTLEIQQQALLREQQKRL NRIKMQEGAKVDLDAIPSAK VREQRMPRDDTSDFLKNSLL ESDSAFIGAYGETYPAIEDD VLPPPSQLPSARERRRNKWK GLDIDSSRPNVAPDGLSLKS ISSVNVDELVRNEERMRL NEFHNKPINTDDESSLVDPD DIMKHIGDDGSNSVATEPWL RPGTSETLKRMAEQLNQE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081