

Anti-GABRB3 antibody (330-460) (STJ112055)

STJ112055

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2562
Gene Symbol	GABRB3
Uniprot ID	GABRB3_HUMAN
Immunogen	
Immunogen Region	330-460
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 330-460 of human GABRB3 (NP_000805.1).
Immunogen Sequence	IFFGRGPGQKQKLAETAKA KNDRSKSESNRVDAHGNIILL TSLEVHNEMNEVSGGIGDTR NSAISFDNSGIQYRKQSMMPR EGHGRFLGDRSLPHKKTHLR RRSSQLKIKIPDLTDVNAID RWSRIVFPFTF



Western blot analysis of extracts of various cell lines, using GABRB3 antibody (STJ112055) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081