

Anti-rprd1b antibody (80-326) (STJ112141)

STJ112141

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:100-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58490
Gene Symbol	RPRD1B
Uniprot ID	RPR1B_HUMAN
Immunogen	
Immunogen Region	80-326
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 80-326 of human rprd1b (Q9NQG5).
Immunogen Sequence	EFESVLVDAFSHVAREADEG CKKPLERLLNIWQERSVYGG EFIQQLKLSMEDSKSPPPKA TEEKSLKRTFQQIQEEEDD DYPGSYSPQDPSAGPLLTEE LIKALQDLENAASGDATVRQ KIASLPQEVQDVSLLKIDT KEAAERLSKTVDEACLLAE YNGRLAAELEDRRQLARMLV EYTNQKQDVLSEKEKKLEFY KQKLARVTQVRKELKSHIQS LPDLSLLPNVTGGLAPLPS



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081