

Anti-RBBP8 antibody (648-897) (STJ112239)

STJ112239

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5932
Gene Symbol	RBBP8
Uniprot ID	CTIP_HUMAN
Immunogen	
Immunogen Region	648-897
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 648-897 of human RBBP8 (NP_002885.1).
Immunogen Sequence	VSFENIQWSIDPGADLSQYK MDVTVIDTKDGSQSKLGGET VDMDCITLVSETVLLKMKKQE QKGEKSSNEERKMNDSLEDML FDRTHTEEYESCLADSFSA ADEEEELSTATKKLHTHGDK QDKVKQKAFVEPYFKGDERE TSLQNFPHIEVVRKKEERRK LLGHTCKECEIYYADMPAEE REKKLASCSRHRFRYIPNT PENFWEVGFSPQTQCMERGY IKEDLDPCPRPKRRQPYNA



Western blot analysis of Rat thymus, using RBBP8 antibody (STJ112239) at 1:1000 dilution. Secondary antibody: HRP-Coat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081