

Anti-CHST6 antibody (196-395) (STJ112267)

STJ112267

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4166
Gene Symbol	CHST6
Uniprot ID	CHST6_HUMAN
Immunogen	
Immunogen Region	196-395
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 196-395 of human CHST6 (NP_067628.1).
Immunogen Sequence	RIVHLVRDPRAVLRSEQTA KALARDNGIVLGTNGTWVEA DPGLRVVREVCRRSHVRIAEA ATKPPPFRLRGYRLVRFED LAREPLAEIRALYFTGLSL TPQLEAWIHNIHGGSGPGAR REAFKTSSRNALNVSQAWRH ALPFAKIRRVQELCAGALQL LGYRPVYSEDEQRNLALDLV LPRGLNGFTWASSTASHPRN



Western blot analysis of extracts of various cell lines, using CHST6 antibody (STJ112267) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081