

Anti-CACNA2D3 antibody (900-1040) (STJ112353)

STJ112353

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55799
Gene Symbol	CACNA2D3
Uniprot ID	CA2D3_HUMAN
Immunogen	
Immunogen Region	900-1040
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1040 of human CACNA2D3 (NP_060868.2).
Immunogen Sequence	LYDYQAMCRANKESSDGAHG LLDPYNAFLSAVKWIMTELV LFLVEFNLCSSWWHSDMTAKA QKLKQTLPCDTEYPAFVSE RTIKETGTNIACEDCSKSFV IQQIPSSNLFMVVDSSCLC ESVAPITMAPIEIRYNESLK C



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081