

Anti-SLC26A6 antibody (509-758) (STJ112361)

STJ112361

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	65010
Gene Symbol	SLC26A6
Uniprot ID	S26A6_HUMAN
Immunogen	
Immunogen Region	509-758
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 509-758 of human SLC26A6 (NP_599025.2).
Immunogen Sequence	QMPHYSLGQVPDIDIYRDV AEYSEAKEVRGVKVFRRSAT VYFANAEFYSDALKQRCGVD VDFLISQKKLLKKQEQLKL KQLQKEEKLKQAASPKGAS VSINVNTSLEDMRSNNVEDC KMMVSSGDKMEDATANGQED SKAPDGSTLKALGLPQPDFH SLILDLGALSFVDTVCLKSL KNIFHDFREIEVEVYMAACH SPVVSQLEAGHFFDASITKK HLFASVHDAVTFALQHPRP



Western blot analysis of extracts of 293T cells, using SLC26A6 antibody (STJ112361) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081