

Anti-LMAN1 antibody (270-480) (STJ112470)
STJ112470

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:100 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	3998
Gene Symbol	LMAN1
Uniprot ID	LMAN1_HUMAN
Immunogen	
Region	270-480
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 270-480 of human LMAN1 (NP_005561.1).
Immunogen Sequence	PGKEPPTPDKEISEKEKEY QEEFEHFQQELDKKKEEFQK GHPDLGGQPAEEIFESVGDR ELRQVFEGQNRHLEIKQLN RQLDMILDEQRRYVSSLTEE ISKRGAGMPGQHGQITQQL DTVVKTQHEILRQVNEMKNS MSETVRLVSGMQHPGSAGGV YETTFHFIDIKHIVKRD IDNLVQRNMPSNEKPKCEP PFPSCSLSTVH




Immunofluorescence analysis of U-2 OS cells using LMAN1 Rabbit polyclonal antibody (STJ112470) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH-3T3 cells using LMAN1 Rabbit polyclonal antibody (STJ112470) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of C6 cells using LMAN1 Rabbit polyclonal antibody (STJ112470) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using LMAN1 Rabbit polyclonal antibody (STJ112470) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081