

Anti-MTA2 antibody (60-180) (STJ112866)

STJ112866

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9219
Gene Symbol	MTA2
Uniprot ID	MTA2_HUMAN
Immunogen	
Immunogen Region	60-180
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human MTA2 (NP_004730.2).
Immunogen Sequence	DSNAREFEEESKQPGVSEQQ RHOLKHRELFLSRQFESLPA THIRGKCSVTLLNETDILSQ YLEKEDCFFYSLVFDPVQKT LLADQGEIRVGCKYQAEIPD RLVEGESDNRNQKMEMKVV D