

Anti-MAP3K14 antibody (330-662) (STJ113190)

STJ113190

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9020
Gene Symbol	MAP3K14
Uniprot ID	M3K14_HUMAN
Immunogen	
Immunogen Region	330-662
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 330-662 of human MAP3K14 (NP_003945.2).
Immunogen Sequence	KFSVEEYLVHALQGSVSSGQ AHSLSLAKTWAARGSRRE PSPKTEDNEGVLLEKLPV DYEYREEVHWATHQLRLGRG SFGVEVRMEDKQTGFQCAVK KVRLEVFRAEELMACAGLTS PRIVPLYGAVREGPWVNIFM ELLEGGSLGQLVKEQGCLPE DRALYYLGQALEGLEYLHSR RILHGDVKADNVLLSSDGS AALCDFGHAVCLQPDGLGKS LLTGDIYIPGTETHMAPEV



Western blot analysis of extracts of Rat lung, using MAP3K14 antibody (STJ113190) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081