

Anti-SEC61A1 antibody (310-420) (STJ113219)
STJ113219

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:100-1:500 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29927
Gene Symbol	SEC61A1
Uniprot ID	S61A1_HUMAN
Immunogen	
Immunogen Region	310-420
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 310-420 of human SEC61A1 (NP_037468.1).
Immunogen Sequence	ARFSGNLLVSLGLTWSDTSS GGPARAYPVGGLCYYLSPPE SFGSVLEDPVHAVVYIVFML GSCAFFSKTWIEVSGSSAKD VAKQLKEQQMVMRGHRETSM VHELNRYIPTA

