

## Anti-CSTF3 antibody (1-103) (STJ113266)

STJ113266

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/ELISA           |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                        |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                                                             |
| <b>Clone ID</b>       |                                                                                                                                                        |
| <b>Concentration</b>  | Lot specific                                                                                                                                           |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                           |
| <b>Purification</b>   | Affinity purification                                                                                                                                  |
| <b>Dilution Range</b> | WB:1:500-1:2000<br>ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>    | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                       |
| <b>Isotype</b>        | IgG                                                                                                                                                    |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                         |
| <b>Instruction</b>    |                                                                                                                                                        |

### TARGET INFORMATION

|                           |                                                                                                                      |
|---------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">1479</a>                                                                                                 |
| <b>Gene Symbol</b>        | <a href="#">CSTF3</a>                                                                                                |
| <b>Uniprot ID</b>         | <a href="#">CSTF3_HUMAN</a>                                                                                          |
| <b>Immunogen</b>          |                                                                                                                      |
| <b>Immunogen Region</b>   | 1-103                                                                                                                |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human CSTF3 (NP_001028677.1). |
| <b>Immunogen Sequence</b> | MSGDGATEQAAEYVPEKVKK AEKKLEENPYDLDAWSILIR EAQNQPIDKARKTYERLVAQ FPSSGRFWKLYIEAEVTILF YFFLYQYCSIHCSDRKQVRN IAN         |



Western blot analysis of extracts of various cell lines, using CSTF3 antibody (STJ113266) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 50s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081