

Anti-LIMA1 antibody (560-759) (STJ113282)

STJ113282

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:100-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51474
Gene Symbol	LIMA1
Uniprot ID	LIMA1_HUMAN
Immunogen	
Immunogen Region	560-759
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 560-759 of human LIMA1 (NP_057441.1).
Immunogen Sequence	DEISKPEVPEDVDLDLKKLR RSSSLKERSRPFTVAASFQS TSVKSPKTVSPPIRKGWSMS EQSEESVGGRAERKQVENA KASKKNGNVGKTTWQNKESK GETGKRSEKHSLEMENENL VENGADSDDDNSFLKQQSP QEPKSLNWSSFVDNTFAEEF TTQNQKSQDVELWEGEVVKE LSVEEQIKRNRYYDEDEDEE



Western blot analysis of lysates from A-549 cells, using LIMA1 Rabbit polyclonal antibody (STJ113282) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ-S000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081