

Anti-HLA-DRA antibody (26-216) (STJ113370)

STJ113370

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	3122
Gene Symbol	HLA-DRA
Uniprot ID	DRA_HUMAN
Immunogen	
Immunogen Region	26-216
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-216 of human HLA-DRA (NP_061984.2).
Immunogen Sequence	IKEEHVIIQAEFYLNPDQSG EFMFDGDEIFHVDMAKKE TVWRLEEFGRFASFEAQGAL ANIADV KANLEIMTKRSNYT PITNVPPETVLTNTPVELR EPNVLCFIDKFTPPVNVNVT WLRNGKPVTTGVSETVFLPR EDHLFRKFHYLPFLPSTEDV YDCRVEHWGLDEPLLKHWEF DAPSPLPETTE




Immunohistochemistry analysis of HLA-DRA in paraffin-embedded rat spleen using HLA-DRA Rabbit polyclonal antibody (STJ113370) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of HLA-DRA in paraffin-embedded rat lung using HLA-DRA Rabbit polyclonal antibody (STJ113370) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of HLA-DRA in paraffin-embedded human liver cancer using HLA-DRA Rabbit polyclonal antibody (STJ113370) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of HLA-DRA in paraffin-embedded human liver using HLA-DRA Rabbit polyclonal antibody (STJ113370) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081