

Anti-NFE2L2 antibody (243-400) (STJ113550)

STJ113550

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000 IHC-P: 1:50-1:100 IF/ICC: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	4780
Gene Symbol	NFE2L2
Uniprot ID	NF2L2_HUMAN
Immunogen	
Immunogen Region	243-400
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 243-400 of human NRF2 (NP_006155.2).
Immunogen Sequence	SPHFLNAFEDSFSSILSTED PNQLTVNSLNSDATVNTDFG DEFYSAFIAEPSISNSMPSP ATLSHSLSELLNGPIDVSDL SLCKAFNQNHPESTAEFNDS DSGISLNTSPSVASPEHSVE SSSYGDITLLGLSDSEVEELD SAPGSVKQNGPKTPVHSS



Immunohistochemistry analysis of NRF2 in paraffin-embedded human placenta using NRF2 Rabbit polyclonal antibody (STJ113550) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of NRF2 in paraffin-embedded human breast cancer using NRF2 Rabbit polyclonal antibody (STJ113550) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of lysates from wild type (WT) and NRF2 knockout (KO) HeLa cells using NRF2 Rabbit polyclonal antibody (STJ113550) at 1:1000 dilution. wild type (WT) and HeLa-NFE2L2 knockout (KO) C6 cells were treated by MG132 (50 µM) at 37 °C for 90 minutes. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 45s.