

Anti-IL1A antibody (113-267) (STJ113557)
STJ113557

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3552
Gene Symbol	IL1A
Uniprot ID	IL1A_HUMAN
Immunogen	
Immunogen Region	113-267
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 113-267 of human IL-1 alpha (NP_000566.3).
Immunogen Sequence	SAPFSFLSNVKYNFMRIIKY EFILNDALNQSIIIRANDQYL TAAALHNLDEAVKFDMGAYK SSKDDAKITVILRISKTLQLY VTAQDEDQPVLLKEMPEIPK TITGSETNLLFFWETHGTKN YFTSAHPNLFIAATKQDYWW CLAGGPPSITDFQIL




Immunohistochemistry analysis of IL-1 alpha in paraffin-embedded human lung cancer using IL-1 alpha Rabbit polyclonal antibody (STJ113557) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of IL-1 alpha in paraffin-embedded human colon carcinoma using IL-1 alpha Rabbit polyclonal antibody (STJ113557) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of recombinant Human IL-1 alpha Protein, using IL-1 alpha Rabbit polyclonal antibody (STJ113557) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 10ng/5ng per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.



Western blot analysis of lysates from Rat lung, using IL-1 alpha Rabbit polyclonal antibody (STJ113557) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081