

Anti-EEF1G antibody (168-437) (STJ113571)

STJ113571

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	1937
Gene Symbol	EEF1G
Uniprot ID	EF1G_HUMAN
Immunogen	
Immunogen Region	168-437
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 168-437 of human eEF1G (NP_001395.1).
Immunogen Sequence	LLWLYKQVLEPSFRQAFNPNTNRWFLTCINQPQFRAVLGEVKLCEKMAQFDAQKFAETQPKKDTPRKEKGSREEKQKPQAE RKEEKKAAPAEPEEMDECEQALAAEPKAKDPFAHLPKSTFVLDEFKRKYSNEDTSLVALPYFWEHFDKDGWSLWYSEYR FPEELTQTFMSONLITGMFQRLDKLRKNAFASVILFGTNNSSISGVVWFRGQELAFPLSPDWQVDYESYTWRLDPPGS




Western blot analysis of extracts of various cell lines, using eEF1G antibody (STJ113571) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunofluorescence analysis of NIH-3T3 cells using eEF1G Polyclonal Antibody (STJ113571) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded rat brain using eEF1G rabbit polyclonal antibody (STJ113571) at dilution of 1:50 (40x lens).



Immunofluorescence analysis of C6 cells using eEF1G Polyclonal Antibody (STJ113571) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081